
















EXTERNAL FIXATION HYBRID

Back to basics,
simplicity, stability
and versatility



IT'S GOT TO MAKE SENSE

COMPONENTS	DESCRIPTION	REFERENCE
	Carbon ring inside diameter 160mm inside diameter 200mm	1020-1002 1020-1003
	Ring clamp	1020-1010
	Kirschner Wire Holder	1020-1020
	Kirschner Wire diameter 1,8mm, length 400mm	1010-1030
	Connection diameter 9mm/9mm	1010-1040
	Reduction sleeve diameter 9mm/5mm	1010-1054
	Carbon rod diameter 9mm, length 70mm diameter 9mm, length 90mm diameter 9mm, length 110mm diameter 9mm, length 150mm diameter 9mm, length 200mm	1010-1062 1010-1063 1010-1064 1010-1065 1010-1066
	Self drilling /self tapping pin Pin diameter 5mm, length 140mm, thread length 50mm Pin diameter 5mm, length 110mm, thread length 40mm	2020-1020 2020-1035
	Screw handle	5015-1000
	Drill brace	5010-1000
	Cutting plyer	5040-1000
	Tensioner	5030-1000
	Sterilisation tray	6010-1005

Distributed by :

CE 0459
EN 93/42



Cat. No: HYB 02-2005



Gexfix S.A.
Av. De La Praille 50
CH-1227 Carouge
Tél : +41 22 308 16 70
Fax : +41 22 308 16 79
www.gexfix.ch
ISO 9001/2000
ISO 13485/2003

To ensure the best quality of its products and their improvements,
Gexfix reserves the right to modify all or part of their products