




# EXTERNAL FIXATION SMALL

Back to basics,  
simplicity, stability  
and versatility



IT'S GOT TO MAKE SENSE

**COMPONENTS**

**DESCRIPTION**

**REFERENCE**




**Connection**  
diameter 9mm/9mm

1010-1040



**Sleeve**  
diameter 9mm/3mm  
diameter 9mm/4mm  
diameter 9mm/5mm

1010-1052  
1010-1053  
1010-1054




**Pin cluster clamp**

1010-1070



**Carbon rod**  
diameter 9mm, length 70mm  
diameter 9mm, length 90mm  
diameter 9mm, length 110mm  
diameter 9mm, length 150mm  
diameter 9mm, length 200mm

1010-1062  
1010-1063  
1010-1064  
1010-1065  
1010-1066



**Lengthening /compression tube**

1010-1100



**Self drilling /self tapping pin**  
Pin diameter 3mm, length 70mm, thread length 20mm  
Pin diameter 4mm, length 100mm, thread length 35mm  
Pin diameter 5mm, length 110mm, thread length 40mm

2010-1020  
2015-1020  
2020-1035



**Trocar and sleeve with handle**

5025-2000



**Screw handle**

5015-1000



**Tray**

6010-1001

Distributed by :

CE 0459  
EN 93/42



Cat. No: SM 02-2005



**Gexfix S.A.**  
Av. De La Praille 50  
CH-1227 Carouge  
Tél : +41 22 308 16 70  
Fax : +41 22 308 16 79  
www.gexfix.ch  
ISO 9001/2000  
ISO 13485/2003

To ensure the best quality of its products and their improvements,  
Gexfix reserves the right to modify all or part of their products