

Caslau[®]

fixation sytem kit

Fast execution time (about 15 min)

Wrist fractures represent approximately 17% of all fractures that enter the Emergency department.

The sterile CASLAU system has as its main indication all the composite or split fractures of the wrist, except meta-epiphyseal comminuted fractures (C3).

The CASLAU system has as its main indication all the composite or split fractures of the wrist, except meta-epiphyseal comminuted fractures (C3).

Advantages over non-operative treatment with plaster

- fast execution time (about 15 min)
- very low risk of infection
- immediate mobilization

Advantages over the tough treatment with external bridge fixator

- immediate mobilization

Advantages compared to treatment with bloody volar plate of dorsal

- fast execution time (about 15 min)
- Perctaneous Synthesis
- very low risk of infection
- immediat mobilization



Product description

- 1 Steel L bar to 5 mm Ø 130 mm x 60 mm
- 2 Kirschner wire, 1.8 Ø x 150 mm
- 2 Self-drilling pin
3 mm Ø x length. 70 mm / thread 35 mm
- 4 connector 9 mm Ø / Ø 9 mm
- 2 Reducer 9 mm / 1.8 mm
- 2 Reducer 9 mm / 3 mm
- 4 Reducer 9 mm / 5 mm


References and variants available 1060-1001

We care
about our clients



Caslau®

fixation system kit



The Gexfix device is a high technology, lightweight, x-ray transparent, external fixation system for trauma orthopaedics. Its design and manufacture make it available to a great many individuals. It is made in Switzerland and distributed worldwide.

The Gexfix fixator affords total freedom to the surgeon to insert the PINS, as well as an optimal ease of use thanks to the reduced number of ancillary items.

Gexfix external fixation orthopaedics

Gexfix SA
Avenue de la Praille 50
CH - 1227 Carouge

T +41(0)22 308 16 70
F +41(0)22 308 16 79

info@gexfix.ch
www.gexfix.ch



Back to basics

simplicity, stability and versatility