



**Gexfix** external fixation orthopaedics

Gexfix SA  
Avenue de la Praille 50  
CH - 1227 Carouge

T +41(0)22 308 16 70  
F +41(0)22 308 16 79

info@gexfix.ch  
www.gexfix.ch



*More than surgery*

# Gexfix

## New products


**We care**  
about our clients



## Gex-Finger

### Description:

The GEX-FINGER kit (REF: 1060-1012) is composed of 2 rods, 2 nuts, 4 springs, 3 K-wires and 6 K-wire caps. The rods and nuts are made of Polyether Ether Ketone (PEEK), a radiolucent, robust thermoplastic biomaterial used in medical implants. The springs and K-wires are made of stainless steel.



## MiniBone

### Description:


Hand and foot fractures.  
Metacarpal bone lengthening and metatarsals.  
Processing brachymétatarsalgies.  
REF: 1060 -1014



## MiniBone facial

### Description:


Facial reconstruction  
Mandibles elongation



## New D.O.S. System Adapted MRI\*

### Description:


The D.O.S. System was created in response to the need for a new minimally invasive reduction and osteosynthesis device for the treatment of fractures of the upper limbs (shoulder, humerus and elbow)  
REF: 1030-1000




## Gex-Foot

### Description:

Treatment of fractures of the foot and ankle



## Easy-clip Connectors



REF: (LARGE) Ø 15:  
1010-1006 (TITANIUM SCREW)  
1010-1005 (S. STELL SCREW)

REF: (SMALL) Ø 9:  
1010-1046 (TITANIUM SCREW)  
1010-1045 (STELL SCREW)

## Titanium Pin & K-Wire ADAPTED FOR MRI\*

\*Not tested